

OI-922

Page 38 of 67 Rev. BS

Vialed and Packaged Test Strip Release

Test Strips

nova
biomedical

CERTIFICATE OF COMPLIANCE

This product has been manufactured in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) standards and meets Nova Biomedical product specifications.

Customer:	N/A	
Description:	TEST STRIPS, STATSTRIP CAPILLARY GLUCOSE ROW, PACKAGED 100 COUNT	
Catalog Number:	61500	
Nova Biomedical P/N	61500	
Target Range (mg/dL)	SD(mg/dL) /CV (%)	
30-50	≤ 4 mg/dL	PASS / FAIL
51-110	≤ 6.0%	PASS / FAIL
111-150	≤ 5.0%	PASS / FAIL
151-250	≤ 5.0%	PASS / FAIL
251-400	≤ 5.0%	PASS / FAIL
Accuracy Criteria 1: Total within +/-12 mg/dL or +/- 12% of reference method ≥ 95%		PASS / FAIL
Accuracy Criteria 2: Total within +/-15 mg/dL or +/- 15% of reference method ≥ 98%		PASS / FAIL
Qty:	3,760	
Lot No.	0326194309	
Mfg. Date:	2026-07-13	
Expiration Date:	2029-01-13	
Unique Device Identifier (UDI):	(01)20385480615000 (11)260713 (17)290113 (10)0326194309	

Mark Browne
Quality Control Manager

7/27/26
Date

Nova Biomedical Corporation

200 Prospect Street Waltham, MA 02454