

Attachment 3



Certificate of Conformance & Analysis

Nova Biomedical certifies that the following Nova Vet Product has been built and tested in compliance with ISO 13485:2016 and Medical Device Single Audit Program (MDSAP) and meets Nova Biomedical product specifications.

Description: Test strip 25ct Novavet Glu
 Part Number: 53069
 BOM Rev.: C
 Lot Number: 1750223088
 Expiration Date: 2025-03-29
 Date of Manufacture: 2023-03-29
 UDI: (01)20385480 530693
 (11) 230329
 (17) 250329
 (10) 1750223088
 Quantity: 301
 Date of Release: 2023-04-14

Release Test Results Summary					
Target Level (mg/dL)	YSI Mean (mg/dL)	Mean (mg/dL)	SD (mg/dL) /CV (%)	SD/CV Specification	P/F
40-80	59.9	66	3.4	≤ 9 mg/dL	P
120-180	167.0	180	2.2	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
220-300	271.5	285	2.0	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
320-440	412.0	413	3.5	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
450-550	490.0	508	1.3	Fail If Both CV of this level and combined CV (L2-L5) >5.0%	P
Total within ISO Zone and %: <u>359/99.7%</u> (per ISO15197:2013)				Total within ISO Zone ≥95%	P

NOTE: Test data in table above taken from completed and approved Data Sheets from TP 139-501 for this lot number.

 Quality Control Manager Signature
 Title: SR Mfg. Manager
 Date: 4/14/23

Nova Biomedical Corporation 200 Prospect Street Waltham, MA 02454