

Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET. Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen. Αυτόματη συσκευή με κρεατινίνη επιπέδου ελέγχου χημείας αέρσιων αίματος, CO-Οξυμέτρο, CO-Οξυμέτρο Stat Profile Prime Plus® VET. Cartușcia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET. Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET con creatinina, Stat Profile Prime Plus® VET vérégáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininnal. (מחשבת ביקורת כימיה ואטומטית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス、CO-オキシメーター、生化学検査用コントロール自動カートリッジ(クレアチニン)、クレアチン入り Stat Profile Prime Plus® VET 혈액 가스, CO-소 농도계, 화학 조제제 자동 카트리지, Stat Profile Prime Plus® VET 気血、一酸化炭素血気、化学对照液自動試劑盒 (含肌酐)

LOT 23214030 2024-12-29

CONTROL 1 2 3 4 5

Table with columns for Analyte, CONTROL 1, CONTROL 2, CONTROL 3, CONTROL 4, and CONTROL 5. Rows include pH, H+, PCO2, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, iMg, Glu, Glu, Lac, Lac, BUN, BUN, Urea, Urea, Creatinine, Creatinine, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, tBil.

EN

ES

FR

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematoct (Hct), fetal hemoglobin (HbF), total hemoglobin (tHb), oxymoglobin (O2Hb), methemoglobin (MetHb), and desoxyhemoglobin (HHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (urea), and Creatinine in Level 4 and 5. A change kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers ONLY.

Descripción del producto: Material acuoso de control de calidad para supervisar el desempeño de pH, PCO2, PO2, SO2, hematocrito (Hct), hemoglobina fetal (HbF), hemoglobina total (tHb), oxihemoglobina (O2Hb), carboxihemoglobina (COHb), metahemoglobina (MetHb) y desoxihemoglobina (HHb) en los niveles 1, 2 y 3, además de Na+, K+, Cl-, Ca, Mg, glucosa, lactato, BUN (urea) y creatinina en los niveles 4 y 5. Se incluye un kit de carga compuesto por dos jeringas llenas de solución de creatinina para activar el cartucho antes de la instalación. Para uso UNICAMENTE con los analizadores Stat Profile Prime Plus VET.

Description du produit: Matériau aqueux de contrôle de qualité destiné à surveiller la performance des éléments suivants : pH, PCO2, PO2, SO2, hématocrite (Hct), hémoglobine totale (tHb), hémoglobine totale (tHb), oxymoglobine (O2Hb), carboxyhémoglobine (COHb), méthémoglobine (MetHb) et désoxyhémoglobine (HHb) dans les niveaux 1, 2 et 3, ainsi que Na+, K+, Cl-, Ca, Mg, glucose et lactate, BUN (urée) et créatine dans les niveaux 4 et 5. Un kit de charge composé de 2 seringues remplies de solution de créatine est inclus pour l'activation de la cartouche avant l'installation. Pour utilisation avec les analyseurs Stat Profile Prime Plus VET UNIQUEMENT.



200 Prospect Street Waltham, MA 02454 U.S.A.

