



REF 54507, 57957, 58102

15°C  30°C

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 23304026

Exp. Date: 2025-04-17

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.343	pH	6.840	PCO ₂	50.9 mmHg	KCL	2M
PCO ₂	25.4 mmHg			PO ₂	99.2 mmHg	 <div>WASTE CONTAINER</div> <div>W</div>	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02