

Stat Profile Prime Plus® Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus®, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® avec créatine, Stat Profile Prime Plus® Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φοισίγγιο υλικού ελέγχου χημείας αέρσιων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus®, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus®, Cartucho automático de controles de química, de CO-oxímetro e de Gás no sangue Stat Profile Prime Plus® com creatinina, Stat Profile Prime Plus® vérgáz, CO-oximéter, kémiai kontrolló auto automatikus patron kreatininnel, קריאטינין, כימיה אוטומטית בקרת תוספית ברת Stat Profile Prime Plus® Co-Oximeter וגז, Stat Profile Prime Plus® 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ (クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® 血気, CO-二酸化炭素血気儀, 化学对照溶液自動試劑盒 (含肌酐)

LOT 24110032

2025-10-09

CONTROL 1 2 3 4 5

Table with 6 columns: Analyte, Control 1, Control 2, Control 3, Control 4, Control 5. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iMg, iGlu, Lac, BUN, Urea, Creatinine, HbF+, tHb, COHb, MetHb, HHb, tBil+, tBil-.

EN Product Description: Adequate quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF+), total hemoglobin (tHb), total bilirubin (tBil), oxymethoglobin (OxHb), carboxyhemoglobin (COHb), methemoglobin (MetHb), and cyanmethemoglobin (CyHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urea), and creatinine in Level 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus Analyzers ONLY.

ES Descripción del producto: Material adecuado de control de calidad para supervisar el desempeño de pH, PCO2, PO2, SO2, hematocrito (Hct), hemoglobina fetal (HbF+), hemoglobina total (tHb), bilirubina total (tBil), oxihemoglobina (OxHb), metahemoglobina (MetHb) y cianohemoglobina (CyHb) en niveles 1, 2 y 3, además de Na+, K+, Cl-, Ca, Mg, glucosa, lactato, BUN (urea) y creatinina en niveles 4 y 5. Se incluye un kit de carga compuesto por dos jeringas llenas de solución de creatinina para activar el cartucho antes de su uso. Para usar UNICAMENTE con los analizadores Stat Profile Prime Plus.

FR Description du produit: Matière de contrôle de la qualité adéquate pour surveiller la performance des éléments suivants: pH, PCO2, PO2, SO2, hématocrite (Hct), hémoglobine totale (tHb), hémoglobine totale (tHb), bilirubine totale (tBil), oxyméthoglobine (OxHb), méthémoglobine (MetHb), cyanométhémoglobine (CyHb), méthémoglobine (MetHb) et cyanométhémoglobine (CyHb) dans les niveaux 1, 2 et 3 ainsi que Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urée) et créatinine dans les niveaux 4 et 5. Un kit de charge composé de 2 seringues remplies de solution de créatinine est inclus pour l'activation de la cartouche avant l'installation. Pour utilisation avec les analyseurs Stat Profile Prime Plus UNIQUEMENT.



