

Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cardridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φασίγγιο με κρεατινίνη επιπέδου ελέγχου χημείας αερίων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetero e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérgáz, CO-oxímeter, kémiai kontrollok automatikus patron kreatininnel, (מחסנית בקרת כימיה אוטומטית עם קראטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cardridge, Stat Profile Prime Plus® VET 血液ガス、CO オキシメーター、生化学検査用コントロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절액 자동 카트리지, Stat Profile Prime Plus® VET 血气、一氧化炭-血氧仪、化学对照溶液自动试剂盒 (含肌酐)

LOT 24117039

CONTROL 1 2 3 4 5

2025-10-17

Table with 7 columns: Parameter, CONTROL 1 (min - X - max), CONTROL 2 (min - X - max), CONTROL 3 (min - X - max), CONTROL 4 (min - X - max), CONTROL 5 (min - X - max), and a final column for expected ranges. Parameters include pH, H+, PCO2, PO2, PO, SO2, Hct, Na+, K+, Cl-, iCa, Ca, iMg, mGlu, Glu, Lac, BUN, Urea, Creatinine, HbF, tHb, tHb, tHb, O2Hb, COHb, MetHb, HHb, tBil, and iBil.

EN

ES

FR

Product Description: Aqueous quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), fetal hemoglobin (HbF), total hemoglobin (tHb), oxyhemoglobin (O2Hb), carboxyhemoglobin (MetHb), and desoxygenhemoglobin (HHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, Glucose, Lactate, BUN (urea), and Creatinine in Levels 4 and 5. A charge kit consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers ONLY.

Descripción del producto: Material acuoso de control de calidad para supervisar el desempeño de pH, PCO2, PO2, SO2, hematocrito (Hct), hemoglobina fetal (HbF), hemoglobina total (tHb), oxihemoglobina (O2Hb), carboxihemoglobina (MetHb) y desoxihemoglobina (HHb) en niveles 1, 2 y 3, además de Na+, K+, Cl-, Ca, Mg, glucosa, lactato, BUN (urea) y creatinina en niveles 4 y 5. Se incluye un kit de carga con dos jeringas con solución de creatinina para activar el cartucho antes de la instalación. Para usar con los analizadores Stat Profile Prime Plus VET.

Description du produit: Matériau de contrôle de qualité aqueux pour surveiller la performance des éléments suivants : pH, PCO2, PO2, SO2, hématocrite (Hct), hémoglobine totale (tHb), hémoglobine totale (HbF), oxyhémoglobine (O2Hb), carboxyhémoglobine (MetHb) et desoxyhémoglobine (HHb) dans les niveaux 1, 2 et 3 ainsi que Na+, K+, Cl-, Ca, Mg, glucose et lactate, BUN (urée) et créatine dans les niveaux 4 et 5. Un kit de charge composé de 2 seringues remplies de solution de créatine est inclus pour l'activation de la cartouche avant l'installation. Pour utilisation avec les analyseurs Stat Profile Prime Plus VET. UNiquement.



200 Prospect Street Waltham, MA 02454 U.S.A.

