



REF 54507, 57957, 58102

15°C  30°C

# BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 24170004

Exp. Date: 2025-12-06

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.339	pH	6.840	PCO <sub>2</sub>	51.3 mmHg	KCL	2M
PCO <sub>2</sub>	26.1 mmHg			PO <sub>2</sub>	97.6 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO<sub>2</sub> and PCO<sub>2</sub>. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.