



REF 54507, 57957, 58102

15°C  30°C

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 24213027

Exp. Date: 2026-01-24

CAL A		700mL	CAL B		220mL	CAL F		>720mL	SOLN R		300mL
pH	7.336		pH	6.840		PCO ₂	51.3 mmHg		KCL	2M	
PCO ₂	26.3 mmHg					PO ₂	99.4 mmHg		 WASTE CONTAINER W		

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02