


REF 54507, 57957, 58102

15°C  30°C

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 24253046

Exp. Date: 2026-02-28

CAL A		700mL		CAL B		220mL		CAL F		>720mL		SOLN R		300mL	
pH	7.346	pH	6.840	PCO ₂	51.6 mmHg	KCL	2M								
PCO ₂	25.7 mmHg			PO ₂	97.9 mmHg	 WASTE CONTAINER W									

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.