

Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemikontrollen, Αυτόματο φουίγγιο με κρεατινίνη επιπέδου ελέγχου de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérgáz, CO-oximéter, kémiai kontrollok automatikus patron kreatininnel, (מומטת ברתם ביחדה אוטומטית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス, CO オキシメーター, 生化学検査用コントロール自動カートリッジ(クレアチニン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-산소 농도계, 화학 조절제 자동 카트리지, Stat Profile Prime Plus® VET 血気, 一酸化炭素血気儀, 化学对照溶液自動試劑盒 (含肌酐)

LOT 24292044

2026-04-01

CONTROL 1 2 3 4 5

Table with columns for parameters (pH, H+, PCO2, etc.) and rows for expected ranges across five control levels. Includes sub-headers in multiple languages: Expected Ranges, Rangos esperados, Plages attendues, Erwartungsbereiche, Αναμενόμενο εύρος, Intervalli previsti, Intervalos previstos, Vártartományok, הצפויים, 予測範囲, 예상 범위, 预期范围值.

EN

Product Description: Alnovus quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), total hemoglobin (tHb), total hemoglobin (tHb), carboxyhemoglobin (COHb), methemoglobin (MetHb), and desoxyhemoglobin (dHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urea), creatinine and creatinine in Level 4 and 5. A change in packaging of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers Only.

ES

Descripción del producto: Material de control de calidad para supervisar el desempeño de pH, PCO2, PO2, SO2, hematócrito (Hct), hemoglobina total (tHb), hemoglobina total (tHb), carboxihemoglobina (COHb), metahemoglobina (MetHb) y desoxihemoglobina (dHb) en los niveles 1, 2 y 3 además de Na+, K+, Cl-, Ca, Mg, glucosa, lactato, BUN (urea), creatinina y creatinina en los niveles 4 y 5. Se incluye un cambio de empaque por dos jeringas llenas de solución de creatinina para activar el cartucho antes de la instalación. Para usar con los analizadores Stat Profile Prime Plus VET.

FR

Description du produit: Matière de contrôle de la qualité destinée pour surveiller la performance des éléments suivants: pH, PCO2, PO2, SO2, hématocrite (Hct), méthéoglobine totale (tHb), méthéoglobine totale (tHb), carboxyhémoglobine (COHb), méthéoglobine (MetHb) et desoxyhémoglobine (dHb) dans les niveaux 1, 2 et 3 ainsi que Na+, K+, Cl-, Ca, Mg, glucose, lactate, BUN (urée), créatinine et créatine dans les niveaux 4 et 5. Un lot de change d'emballage de 2 seringues remplies de solution de créatine est inclus pour l'activation de la cartouche avant l'installation. Pour utilisation avec les analyseurs Stat Profile Prime Plus VET. Utilisation prévue: Destinée à la surveillance de la performance des analyseurs Stat Profile Prime Plus VET.



200 Prospect Street Waltham, MA 02454 U.S.A.

