


BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 24309089

Exp. Date: 2026-04-23

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.349	pH	6.840	PCO ₂	51.2 mmHg	KCL	2M
PCO ₂	24.7 mmHg			PO ₂	102.2 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.