


REF 54507, 57957, 58102

15°C - 30°C

# BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 24339080

Exp. Date: 2026-05-14

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.350	pH	6.840	PCO <sub>2</sub>	50.7 mmHg	KCL	2M
PCO <sub>2</sub>	24.3 mmHg			PO <sub>2</sub>	98.7 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO<sub>2</sub> and PCO<sub>2</sub>. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.