


REF 54507, 57957, 58102

15°C - 30°C

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 25017068

Exp. Date: 2026-07-06

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.347	pH	6.840	PCO ₂	51.0 mmHg	KCL	2M
PCO ₂	24.5 mmHg			PO ₂	100.7 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.