

Stat Profile Prime Plus® VET Blood Gas, CO-Oximeter, Chemistry Controls Auto-Cartridge with Creatinine

Cartucho automático con creatinina para controles de química, gases en sangre y cooxímetro Stat Profile Prime Plus® VET, Cartouche automatique de contrôles de gaz du sang/CO-oxymètre, chimie Stat Profile Prime Plus® VET avec créatine, Stat Profile Prime Plus® VET Auto-Kassette mit Kreatinin für Blutgas-, CO-Oximeter- und Blutchemiekontrollen, Αυτόματο φύσιγιο με κρεατινίνη επιπέδου ελέγχου χημείας, αερίων αίματος, CO-Οξυμέτρο Stat Profile Prime Plus® VET, Cartuccia con creatinina per controlli automatici chimici per gas ematici/CO-ossimetria Stat Profile Prime Plus® VET, Cartucho automático de controles de química, de CO-oxímetro e de gás no sangue Stat Profile Prime Plus® VET com creatinina, Stat Profile Prime Plus® VET vérzág, CO-oximéter, kémiai kontrolok automatikus patron kreatininnel, (מוסנית בקרת כימיה אוטומטית עם קריאטינין) Stat Profile Prime Plus® VET Chemistry Controls Auto-Cartridge, Stat Profile Prime Plus® VET 血液ガス、CO オキシメーター、生化学検査用コントロール自動カートリッジ (クリアナチン), 크레아티닌 사용 Stat Profile Prime Plus® VET 혈액 가스, CO-오스 도측계, 화학 조절용 자동 카트리지가, Stat Profile Prime Plus® VET 血気、一산화탄소血気仪、化学对照液用自动试剂盒 (舍副劑)

LOT 25042035

CONTROL 1 2 3 4 5

2026-07-24

Table with columns for Analyte, Unit, and Control Ranges (min - X - max) for 5 control levels. Rows include pH, H+, PCO2, PO2, SO2, Hct, Na+, K+, Cl-, iCa, iCa, iMg, Glu, Lac, BUN, Urea, Creatinine, HbF, Hb, iHb, COHb, Methb, HbH, iBil, and Bil.

EN

Product Description: Accurately quality control material for monitoring the performance of pH, PCO2, PO2, SO2, hematocrit (Hct), total hemoglobin (tHb), total hemoglobin (tHb), oxygen hemoglobin (O2Hb), methemoglobin (MetHb), and deoxyhemoglobin (dHb) in Levels 1, 2 and 3 as well as Na+, K+, Cl-, iCa, Mg, Glucose, Lactate, BUN (Urea), and Creatinine in Level 4 and 5. A cartridge consisting of 2 syringes filled with creatinine solution is included for cartridge activation prior to installation. For use with Stat Profile Prime Plus VET Analyzers ONLY.

ES

Descripción del producto: Material de control de calidad para supervisar el desempeño de: pH, PCO2, PO2, SO2, hematocrito (Hct), hemoglobina total (tHb), hemoglobina total (tHb), oxihemoglobina (O2Hb), metahemoglobina (MetHb) y desoxihemoglobina (dHb) en niveles 1, 2 y 3, además de Na+, K+, Cl-, iCa, Mg, glucosa, lactato, BUN (urea) y creatinina en niveles 4 y 5. Se incluye un cartucho con 2 jeringas con solución de creatinina para activar el cartucho antes de la instalación. Para usar ÚNICAMENTE con los analizadores Stat Profile Prime Plus VET.

FR

Description du produit: Matière de contrôle de la qualité aqueux pour surveiller la performance des éléments suivants : pH, PCO2, PO2, SO2, hématocrite (Hct), héoglobine totale (tHb), héoglobine totale (tHb), oxihémoglobine (O2Hb), méthémoglobine (MetHb) et désoxyhéoglobine (dHb) dans les niveaux 1, 2 et 3 ainsi que Na+, K+, Cl-, iCa, Mg, glucose, lactate, BUN (urée) et créatinine dans les niveaux 4 et 5. Un kit de charge composé de 2 seringues remplies de solution de créatine est inclus pour l'activation de la cartouche avant installation. Pour utilisation avec les analyseurs Stat Profile Prime Plus VET UNIQUEMENT.



200 Prospect Street Waltham, MA 02454 U.S.A.

