


REF 54507, 57957, 58102

15°C  30°C

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 25050067

Exp. Date: 2026-08-03

CAL	A	700mL	CAL	B	220mL	CAL	F	>720mL	SOLN	R	300mL
pH	7.346		pH	6.840		PCO ₂	50.9	mmHg	KCL	2M	
PCO ₂	24.6	mmHg				PO ₂	100.5	mmHg		WASTE CONTAINER	W

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.