



CERTIFICATE OF ANALYSIS

BioProfile FAST CDV Reagent Cartridge
Lot Number 25063101

Ref Number 64739
Exp Date 2025-11-15
Date of Manufacture 2024-11-15

CDV Reagent Cartridge L/N:	25063099	Exp Date:	2025-11-15
SOLN SF L/N	24324052	SOLN DP L/N	24320035
	<u>Spec</u>		<u>Result</u>
pH	6.650 - 7.350	Depro Soln	N/A
	Pass		Present
SOLN C L/N	24320034		
	<u>Spec</u>		<u>Result</u>
3% Detergent	N/A		Present

CDV Bottle Pack L/N:	25063100	Exp Date:	2025-11-19
SOLN T L/N	24324032		
	<u>Spec</u>		<u>Result</u>
0.4% Trypan	N/A		Present

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A in with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Haneen Aldarwish

Approval

QC Inspector
Title

3/7/2025
Date