

**BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge**

Lot: 25077062

Exp. Date: 2026-09-07

|                  |           |             |       |                  |           |                                                                                                      |    |
|------------------|-----------|-------------|-------|------------------|-----------|------------------------------------------------------------------------------------------------------|----|
| CAL A 700mL      |           | CAL B 220mL |       | CAL F >720mL     |           | SOLN R 300mL                                                                                         |    |
| pH               | 7.346     | pH          | 6.840 | PCO <sub>2</sub> | 51.9 mmHg | KCL                                                                                                  | 2M |
| PCO <sub>2</sub> | 24.6 mmHg |             |       | PO <sub>2</sub>  | 97.9 mmHg |  WASTE CONTAINER W |    |
|                  |           |             |       |                  |           |                                                                                                      |    |

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO<sub>2</sub> and PCO<sub>2</sub>. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.