

BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge**Lot:** 25108028**Exp. Date:** 2026-09-30

CAL A	700mL	CAL B	220mL	CAL F	>720mL	SOLN R	300mL
pH	7.334	pH	6.840	PCO ₂	53.2 mmHg	KCL	2M
PCO ₂	25.3 mmHg			PO ₂	99.2 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02