


BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 25175029

Exp. Date: 2026-11-19

CAL A 700mL	CAL B 220mL	CAL F >720mL	SOLN R 300mL
pH 7.349 <i>PCO₂</i> 24.2 mmHg	pH 6.840	<i>PCO₂</i> 53.2 mmHg <i>PO₂</i> 97.5 mmHg	KCL 2M
			 WASTE CONTAINER W

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, *PO₂* and *PCO₂*. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02