


BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge

Lot: 25238026

Exp. Date: 2027-01-15

| | | | | | | | |
|--------------------|-----------|--------------------|-------|---------------------|-----------|--|----|
| CAL A 700mL | | CAL B 220mL | | CAL F >720mL | | SOLN R 300mL | |
| pH | 7.332 | pH | 6.840 | PCO ₂ | 52.8 mmHg | KCL | 2M |
| PCO ₂ | 25.1 mmHg | | | PO ₂ | 99.2 mmHg |  WASTE CONTAINER W | |
| | | | | | | | |

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO₂ and PCO₂. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02