


**BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge**

Lot: 26071023

Exp. Date: 2027-08-24

CAL A 700mL		CAL B 220mL		CAL F >720mL		SOLN R 300mL	
pH	7.339	pH	6.840	PCO <sub>2</sub>	53.7 mmHg	KCL	2M
PCO <sub>2</sub>	24.9 mmHg			PO <sub>2</sub>	100.2 mmHg	 WASTE CONTAINER W	

The BioProfile® FLEX 2™ pH/Gas Calibrator Cartridge is for the quantitative determination of pH, PO<sub>2</sub> and PCO<sub>2</sub>. Store at 15°C - 30°C. The waste bag contains preservative, which is an eye, skin and respiratory irritant. Consult the reference manual for directions for use.



200 Prospect Street  
Waltham, MA 02454 U.S.A

LPN 58226B 2020-02