



**CERTIFICATE OF ANALYSIS**

BioProfile FLEX 2 CDV Reagent Pack  
Lot Number 26096050

Ref Number 55297  
Exp Date 2027-02-27  
Date of Manufacture 2026-02-27

	CDV Bottle Pack L/N:	26096049	Exp Date:	2027-02-27	
Soln SF L/N	26064021		Soln DP L/N	26064023	
	<u>Spec</u>	<u>Result</u>		<u>Spec</u>	<u>Result</u>
pH	6.900 - 7.100	Pass	Depro Soln	N/A	Present
Cal C L/N	26064022		Soln T L/N	26058039	
	<u>Spec</u>	<u>Result</u>		<u>Spec</u>	<u>Result</u>
3% Detergent	N/A	Present	0.4% Trypan	N/A	Present

CDV Reagent Cartridge L/N: 26096047      Exp Date: 2027-06-12  
 Soln PF L/N      26071017  
                          Spec      Pass/Fail  
 pH      6.900 - 7.100      Pass

This certifies that this product was manufactured and tested at Nova Biomedical Corporation, Waltham MA 02454 U.S.A in with ISO 13485:2016 Medical Devices Quality Management Systems Requirements, Medical Device Single Audit Program and conforms to the indicated test specifications. The following analytes are traceable to NIST SRM Materials: pH.

Note: Acceptance specifications for this part number are lot dependent and subject to change by the manufacturer

Approval

QC Inspector  
Title

4/15/2026  
Date